



Next-level discovery

Setting a new standard for short reads with the Onso[™] system



15× higher accuracy than other benchtop sequencers



Unprecedented sensitivity



Less sequencing reduces cost + allows higher throughput



Seamless workflow integration

Platform features

Powered by the innovative **sequencing by binding (SBB®) technology**, the Onso system delivers unprecedented accuracy and sensitivity to advance your critical and groundbreaking research.



Higher accuracy

Unique SBB technology delivers 15× higher accuracy for sensitive and specific characterization of low complexity regions with 90% Q40+. Fewer errors mean less false positives and more biological insight.



Streamlined workflow

Simple workflows easily integrate into existing laboratory processes, reducing time and labor.

Improved sensitivity


Using fewer reads, the Onso system lowers the limit of detection, which drives the sensitivity of applications like liquid biopsy. This enables the identification of variants missed by other shortread sequencing platforms.



Efficiency

Q40+ delivers high quality results at lower read depths, enabling scalable sequencing while reducing time and cost. Reduced need for unique molecular identifiers (UMIs) lowers workflow complexity and cost.

Industry-leading sensitivity for low-frequency allele detection in liquid biopsy*



* Analytical validation (TP53 plasmid spike-in) of hybrid capture-based NGS genomic profiling of cell-free circulating tumor DNA † Variant allele frequency (VAF) limit of detection (LoD) at 99.99% specificity for >90% of the sites

The Onso system

The Onso system is a scalable and flexible benchtop platform that gives you the flexibility to integrate existing short-read tools.



Onso system key accessories:

- Library conversion kit
- PacBio® native library prep kits (HMW and fragmented DNA preparation)
- Blocking oligos for hybridization capture methods
- Library amplification kit
- Library quantification kit
- Lambda sequencing control library (spike-in library)

Onso reagents	Read length	Reads	Output (Gb)	Run time	Quality score
200 cycle sequencing kit	2 × 100 bp 1 × 200 bp	400-500 M	100	32 hours	≥90% Q40
300 cycle sequencing kit	2 × 150 bp	400-500 M	150	48 hours	

What can you do with the Onso system?



Cancer research

Enable development of screening, monitoring, and therapy selection



Gene editing

Confirmation of editing outcomes and biomarker discovery



Exome/panels

Validated WES using Twist Bioscience hybrid capture protocol



Single cell

10x Genomics Chromium Next GEM Single Cell 3' libraries supported





Headquarters 1305 O'Brien Drive Menlo Park, CA 94025 United States Phone: 1.650.521.8000 Customer service 1.877.920.PACB (7222), option 1 Fax: 1.650.618.2699 orders@pacb.com Technical support 1.877.920.PACB (7222), option 2 support@pacb.com

As a global company, we have offices in countries around the world. Visit **pacb.com/contact-us** for contact info.

READY TO GET STARTED WITH THE ONSO SYSTEM?

Products and services pacb.com/products

pacb.com/pubs

Documentation pacb.com/documentation Application-specific workflows pacb.com/applications Publications using SMRT® sequencing Connect with PacBio for more info: North America: nasales@pacb.com South America: sasales@pacb.com EMEA: emea@pacb.com Asia Pacific: apsales@pacb.com

Contact a certified service provider **pacb.com/CSP**

Distributed by





Research use only. Not for use in diagnostic procedures. © 2022 Pacific Biosciences of California, Inc. ("PacBio"). All rights reserved. Information in this document is subject to change without notice. PacBio assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of PacBio products and/or third-party products. Refer to the applicable PacBio terms and conditions of sale and to the applicable terms at pacb.com/license. Pacific Biosciences, the PacBio logo, PacBio, Circulomics, Omniome, SMRT, SMRTbell, Iso-Seq, Sequel, Nanobind, and SBB are trademarks of PacBio. All other trademarks are the sole property of their respective owners.

